

IMPORTANT Information



- Do not squeeze, stress or burn a tick – it may regurgitate or inject more saliva into you.
- Do not use petroleum jelly, oil or any chemical products to remove it.
- Never pull or twist a tick as its mouthparts may break off and remain in the skin.

These days, more and more ticks are infected, so it is important to remove ticks safely, quickly and correctly!



Ticks can carry dangerous substances in their saliva and body fluids. Avoid touching removed ticks. Dispose of ticks safely. Do not squash them between your finger nails. If possible disinfect the bite area after removal and wash your hands and the TickCard with hot soapy water.

For more information on ticks and tick diseases see: www.TickCard.co.uk

We support Lyme Disease Charities.
The TickCard and all packaging is designed with quality, durability and the environment in mind. Plastics are 100% PVC-free, with low environmental impact (ABS and PP). Retain this card for reference.



© copyright 2011 PANOSUN Ltd.



made in the EU
Patented: EP1424917

Keep your TickCard™ ready in your wallet



www.TickCard.co.uk
All rights reserved.


10 million sold worldwide

NEW!
now with a magnifying lens

Recommended by Experts


TickCard™

Protect NOW from TICK DISEASE – don't delay proper removal!



Remove ticks safely, quickly & correctly

With a TickCard in your wallet you are well prepared.



credit card size



“Absolutely PERFECT!”
orn Lager, DK, Tick Removal Specialist

Developed for your safety

Don't delay proper removal:

- To minimise the risk of dangerous infections remove ticks safely, correctly and as soon as possible!
- Check yourself and your family for ticks after outdoor pursuits.
- Keep your TickCard handy – don't delay proper removal.

Protect yourself – stay safe!



Keep your TickCard available

Conventional methods of removing ticks can be dangerous.

Biologists found that squeezing, twisting or in any way stressing a tick increases the risk of it injecting harmful agents into the bloodstream of its host, thus potentially causing serious disease. Developed in co-operation with scientists and doctors specialising in tick-borne diseases, the patented TickCard utilises an entirely new principle: ticks are removed safely, quickly and without the use of force. This method helps avoid ticks releasing their dangerous microbes into the bloodstream.

INSTRUCTIONS:

1. Slip the notch of the TickCard under the tick.
2. Gently push the card forward and upward with steady even pressure.
3. Allow the tick to let go and be pulled out of the skin easily and effortlessly.

Remove small, not yet engorged ticks and bee-stings with the smaller notch – it can also be used for thorns and splinters.



Serious tick-borne diseases include *Lyme borreliosis* which is increasingly common throughout Europe (incl. UK) and America. Symptom patterns vary from person to person. Early indications in the days after the bite can include flu-like symptoms, feeling unwell, a rash, headache, muscle pain, tender glands and sensitivity to temperature, sound and light.

Please seek medical advice promptly if you are concerned or experience any of the above symptoms!

Please note: a small red itchy patch is a normal body-reaction to a tick bite.

With handy Magnifying Lens

TickCard™


Also for bee-stings, thorns and splinters

A TickCard can help you avoid potentially serious disease!



Removes big and small ticks – safe and easy to use.

Designed for use anywhere on the body.



With a TickCard™ you are always well prepared.